

SPECIFICATIONS

ROMED CLEANSING WIPES

REF: ABW



Identification

Product name: ROMED cleansing wipes
Ref: ABW

Materials

| | |
|---------------------|-----------------------------------|
| Fabric | Non-woven (thermbond) impregnated |
| Weight | 53 g/m ² |
| Dry weight | 18 g/m ² |
| Dimensions | 200 ± 10 mm x 110 mm ± 10 mm |
| Active ingredient | 80% ethanol, 5% isopropanol |
| Inactive ingredient | Water |
| pH | 6.5 – 7.5 |

Purpose

For the cleaning of small surfaces (~0.5 m²).

Mode of application

Take out a single tissue and wipe the objects to be cleaned. No additional wiping necessary.

Caution

- The tissues are pre-moisturized with an alcoholic liquid which is easily inflammable.
- Ensure that the container is always closed tightly after use.
- Be careful with surfaces sensitive to alcohol (e.g. Plexiglas).
- Avoid contact with eyes. If contacted, flush eyes with copious quantities of water.
- If taken in, DO NOT induce vomiting, but obtain medical attention immediately.
- Stop use if irritation and redness develop on the skin. Ask a doctor if conditions persist for more than 72 hours.

Never use wipe twice

Flammable, keep away from fire or flame

For external use only

Keep out of reach of children

Do not use after expiry date

Single use, discard after use

Keep at room temperature

Packing

110 wipes in a dispenser, 6 dispensers in an inner box, 6 inner boxes in an export carton